


AEROpart®

Oxygen-accessories



AEROpart®

Oxygen-accessories

AEROpart® - Oxygen masks | Inhalation sets

The soft material (vinyl) guarantees a maximum comfort. Three or six channel kink safe tubing and connection adapters assure high security during the application. The strap or ear clip of the nasal cannulas as well as the adjustable nose clip of the masks ensure a comfortable fit. The aerosol inhalation sets are additionally equipped with a standard nebulizer.

Options

For children and adults, with nose clip
and 2 m / 2.1 m kink safe tubing / transparent/clear

O₂-Masks

Without Reservoir



With O₂-Reservoir




Aerosol-Inhalation sets



AEROp^{art}[®] - O₂-Cannulas

Premature infant | Neonates | Infants | Toddler | Children | Adults

Also available with 0.33 m O₂- safety tubing for a patient-side application with FireSafe™ - adaptors




ideal for long-term use

Options

Straight nasal tip Soft

Transparent/clear



with 2 m safety tubing




With 1.5 m safety tubing
incl. ear clip

Curved nasal tips Supersoft

Transparent/clear



with 2 m safety tubing



AEROp^{art}[®] - O₂-Tubes | Adapters / Connectors

Safety tubing, clear, available in different lengths

Options

O₂-Safety tubing

Safety tubing, clear

 available in following lengths 0.33 / 2.1 / 7.5 / 10 / 15 m

 available in following lengths 0.33 / 2.0 / 3.0 / 7.5 / 10 / 15 m



  3 star lumen



 6 star lumen

Adapters/Connectors



Barbed hose adapter



Angled barbed hose adapter



Firesafe™
tubing connector



T-Shape tubing connector



Y-shape tubing connector



Tubing connector
360° rotatable



tubing connector fix



Firesafe™
barbed hose adapter



Water trap

AEROpart®

Oxygen-accessories

AEROpart® - Humidifiers

Humidifiers for use in oxygen therapy.

Different models:

Different capacities | Inlet UNF 9/16" or G 3/8" | Disposable or reusable



Disposable 140 ml
UNF 9/16"




Disposable 350 ml
UNF 9/16"





Autoclavable 250 ml
UNF 9/16"
(also in G 3/8" available)

AEROpart® - Sterile water

Sterile, pyretogeneticfree water, 1000 ml for Inhalation and humidification.



Refill bottle
for manual refilling



With sterile lid for usage
with auto feed systems



AEROp[®] - Holder for Sterile Water Systems and Humidifiers




Holder of transparent acrylic glass for humidifiers or sterile water systems.

Consisting of basic box and exchangeable template covers | With improved platform of 13 x 13 cm | Additional eyelets for fixing

Available either as set consisting of basic box and template cover or individually.


Magnetic application



Basic box



Template cover for AEROp[®] - AB1-916E
(humidifier 140 ml)




Example




HSTW01-H-916E

Template cover for AEROp[®] - AB1-916
and AB1-916HF (humidifier 350 ml)



Aquapak 340

Template cover for sterile water systems Respiflo 325 ml,
Aquapak 340 and 650 ml, STER-O₂ 350 and 550 ml



AEROpart®

Oxygen-accessories

Product overview - Oxygen accessories

AEROpart® - Oxygen masks and Aerosol Inhalation Sets - transparent, clear

AEROpart® - SM01 Oxygen masks with 2.1 m safety tubing, phthalat-free

Item-No.	Description
HSM01-E	O ₂ -Mask for adults
HSM01-K	O ₂ -Mask for children
HSM01-EB	O ₂ -Mask for adults including O ₂ -Reservoir



AEROpart® - AI11 Aerosol Inhalation Sets with mask, nebulizer and 2 m safety tubing

Item-No.	Description
HAI11-E	Aerosol inhalation set for adults
HAI11-K	Aerosol inhalation set for children

AEROpart® - SM11 Oxygen masks with/without 2.1 m safety tubing, with/without reservoir

Item-No.	Description
HSM11-E	O ₂ -Mask for adults
HSM11-K	O ₂ -Mask for children
HSM11-EB	O ₂ -Mask for adults with O ₂ -Reservoir
HSM11-KB	O ₂ -Mask for children with O ₂ -Reservoir
HSM11-EOS	O ₂ -Mask for adults without tube
HSM11-KOS	O ₂ -Mask for children without tube



AEROpart® - O₂-Cannulas - transparent, clear

AEROpart® - SB01 O₂-Cannulas, straight nasal tip, soft (incl. ear clip or phthalat-free)

Item-No.	Description
HSB01-E	O ₂ -Cannulas for adults with 1.5 m safety tubing incl. ear clip
HSB01-S	O ₂ -Cannulas for adults with 2.1 m safety tubing, phthalat-free



AEROpart® - SB02 O₂-Cannulas, curved nasal tip, phthalat-free, supersoft

Item-No.	Description
HSB02-S-0.33	O ₂ -Cannula for adults with 0.33 m safety tubing
HSB02-S	O ₂ -Cannula for adults with 2.1 m safety tubing
HSB02-S-K	O ₂ -Cannula for children with 2.1 m safety tubing
HSB02-S-KK	O ₂ -Cannula for infants with 2.1 m safety tubing
HSB02-S-NS	O ₂ -Cannula for neonates/toddler with 2.1 m safety tubing
HSB02-S-LF	O ₂ -Cannula for adults with 2.1 m safety tubing, for low flows to 3 l/min
HSB02-S-Q	O ₂ -Cannula for adults with 2.1 m safety tubing, for mittlere Flows to 8 l/min
HSB02-S-HF	O ₂ -Cannula for adults with 2.1 m safety tubing, for flows 6 to 8 l/min



AEROpart® - SB11 O₂-Cannulas, straight nasal tip phthalat-free, soft

Item-No.	Description
HSB11-S-0.33	O ₂ -Cannula for adults with 0.33 m safety tubing
HSB11-S	O ₂ -Cannula for adults with 2 m safety tubing
HSB11-S5	O ₂ -Cannula for adults with 5 m safety tubing
HSB11-S-K	O ₂ -Cannula for children with 2 m safety tubing
HSB11-S-KK	O ₂ -Cannula for infants with 2 m safety tubing



AEROpart® - SB12 O₂-Cannulas, curved nasal tip, soft

Item-No.	Description
HSB12-S	O ₂ -Cannula for adults with 2 m safety tubing
HSB12-S-K	O ₂ -Cannula for children with 2 m safety tubing
HSB12-S-KK	O ₂ -Cannula for infants with 2 m safety tubing
HSB12-S-NS	O ₂ -Cannula for neonates/toddler with 2 m safety tubing
HSB02-EZW	Earpads for O ₂ -Cannulas



AEROpart® - O₂-Safety tubings - transparent, clear, with connection pieces

AEROpart® - SS01 O₂-Safety tubings, phthalat-free

Item-No.	Description
HSS01-0.33 / 2.1 / 7.5 / 10 / 15	O ₂ -Safety tubings, different lengths from 0.33 to 15 m



AEROpart® - SS01 O₂-Safety tubings

Item-No.	Description
HSS11-0.33 / 2 / 3 / 7.5 / 10 / 15	O ₂ -Safety tubings, different lengths from 0.33 to 15 m



AEROpart® - Barbed hose adapter, -connector for 6 mm safety tubings

Item-No.	Description
HSA02-6	Barbed hose adapter, white, disposable, UNF 9/16"
HSA02-6G	Barbed hose adapter, green, disposable, UNF 9/16"
HSA02-6-WA	Angled barbed hose adapter, white, disposable, UNF 9/16" outer screw thread
HSA03-6-FS	Firesafe™ barbed hose adapter, UNF 9/16"
HSV01-6	Tubing connector, straight
HSV01-6D	Tubing connector, straight, 360° rotatable
HSV01-6-T	T-Shape tubing connector
HSV01-6-Y	Y-Shape tubing connector
HSV03-6-FS	Firesafe™ tubing connector



Product overview - Oxygen accessories

AEROp[®] - Humidifiers | Sterile water | Adapters / Connectors

AEROp[®] - AB Humidifiers, disposable or autoclavable

Item-No.	Description
HAB01-916	Humidifier, 350 ml, disposable, UNF9/16"
HAB02-916	Humidifier, 250 ml, autoclavable, UNF9/16"
HAB02-38	Humidifier, 250 ml, autoclavable, G 3/8"
HAB01-916-E	Humidifier 140 ml, disposable, UNF9/16"



AEROp[®] - STW 1000 Sterile water

Item-No.	Description
HSTW03-1000FLD	Sterile water, refill bottle for manual refilling, 1000 ml
HSTW03-1000FLDG	Sterile water, with sterile lid, 1000 ml
HWF01	Water trap for oxygen therapy, connectors for 6 mm safety tubings



AEROp[®] - Holder (Basic Box and Template Cover)

Item-No.	Suitable for
HSTW01-H-916E	AEROp [®] - AB1-916E (Humidifier 140 ml)
HSTW01-H-916	AEROp [®] - AB1-916 und AB1-916HF (Humidifier 350 ml)
HSTW01-H-GSU1	Sterile water systems RespiFlo 325 ml, Aquapak 340 and 650 ml, STER-O ₂ 350 and 550 ml



AEROp[®] - Halterungen (Grundbox und Schablone)

Item-No.	Description
HSTW01-H-S916E	Template cover for AEROp [®] - AB1-916E (humidifier 140 ml)
HSTW01-H-S916	Template cover for AEROp [®] - AB1-916 and AB1-916HF (humidifier 350 ml)
HSTW01-H-SGSU1	Template cover for sterile water systems RespiFlo 325 ml, Aquapak 340 and 650 ml, STER-O ₂ 350 and 550 ml
HSTW01-H-BOX	Basic box



HUM

HUM Gesellschaft für Homecare
und Medizintechnik mbH

Zum Pier 79 | D-44536 Lünen

FON +49 (0)2 31/ 88 08 85-0 | FAX +49 (0)2 31/ 88 08 85-58

Mail: sales@hum-online.de | <http://www.hum-online.de>

AEROp[®] - Oxygen-accessories

PP-AEROPART-SZ_en | Status 19.10.2016 | Revision 9

Changes and errors reserved. Illustrations are approximate.

© HUM Gesellschaft für Homecare und Medizintechnik mbH