



ResMed

**Actual size: 140mm (Width)
x 175mm (Height)**

Ruler is for scaling purposes.

**Use this fitting template as
a guide to size your patients.**

How to use the fitting template:


- Cut out the template along the dotted line
- Follow directions on template to determine the recommended cushion size.
- Please ensure the template is printed to scale.

Quattro™ Air

AirFit™ F10

Full face mask
Full face mask **for Her**


Rx Only 
0123




Quattro™ Air


AirFit™ F10

Full face mask
Full face mask **for Her**







SIZE



Sample size check



Rx Only 
0123

