

$R_{\!\!X\,\text{Only}}$

Mask/Device Compatibility List

English

This document specifies the masks compatible with ResMed's positive airway pressure devices and the mask setting, if one is required. For full details regarding the correct use of these devices, please refer to the specific User Guides.

Note: Not all products are available in all regions.

Full Face Masks

Device	AcuCare [™] F1-4	AirFit [™] F30i ¹	AirFit [™] F30	AirFit [™] F20 / AirFit [™] F20 for Her AirTouch [™] F20 / AirTouch [™] F20 for Her	AirFit [™] F10 / AirFit [™] F10 for Her	Mirage Quattro™	Quattro [™] FX / Quattro [™] FX for Her Quattro [™] Air / Quattro [™] Air for Her	Mirage Liberty [™]	Ultra Mirage [™]	Hospital Full Face
AirSense [™] 11	Χ	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face
AirMini™	Χ	Χ	✓	~	Χ	Χ	X	Χ	Χ	Χ
AirStart [™] 10 / AirSense [™] 10	N/V	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face
AirCurve [™] 10 VAuto / AirCurve [™] 10 S / AirCurve [™] 10 ST	N/V	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face
AirCurve [™] 10 ASV / AirCurve [™] 10 ASVAuto	N/V	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face
AirCurve [™] 10 ST-A	N/V	N/V	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face
S9 AutoSet [™] / S9 AutoSet [™] for Her / S9 Elite [™]	N/V	Full Face	N/V	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face
S9 Escape [™] / S9 Escape [™] Auto	N/V	Full Face	N/V	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face
C-Series Tango™	N/V	N/V	N/V	N/V	No	setting requi	red. Consult your clinician before changing	to a different m	nask type.	
S7 Elite [™]	N/V	N/V	N/V	N/V	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL
S8 AutoSet Vantage [™]	N/V	N/V	N/V	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL
S8 Compact [™]	N/V	N/V	N/V	N/V	No	setting requi	red. Consult your clinician before changing	to a different m	nask type.	
S8 Elite [™] / S8 Elite [™] II	N/V	N/V	N/V	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL
S8 Escape [™]	N/V	N/V	N/V	N/V	Requir	es pressure ca	libration at 12 cm H2O. Consult your clinicia	an and check th	ne user guide.	
S8 Escape [™] II	N/V	N/V	N/V	Requi	es pressure calibration a	it 12 cm H2O.	Consult your clinician and check the user gu	uide.		
S8 AutoSet [™] II	N/V	N/V	N/V	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL
VPAP™ III	N/V	N/V	N/V	N/V	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL
VPAP™ III ST-A	N/V	N/V	N/V	N/V	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL
VPAP [™] Adapt SV	N/V	N/V	N/V	N/V	FULL FACE	FULL FACE	FULL FACE	FULL FACE ¹	FULL FACE	FULL FACE
VPAP [™] Adapt (S9)	N/V	Full Face	N/V	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face
VPAP™ Adapt (S9 ASV)	N/V	Full Face	N/V	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face
VPAP™ Auto 25	N/V	N/V	N/V	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL
VPAP [™] Auto (S9)	N/V	Full Face	N/V	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face

Device	AcuCare [™] F1-4	AirFit [™] F30i ¹	AirFit [™] F30	AirFit [™] F20 / AirFit [™] F20 for Her AirTouch [™] F20 / AirTouch [™] F20 for Her	AirFit [™] F10 / AirFit [™] F10 for Her	Mirage Quattro [™]	Quattro [™] FX / Quattro [™] FX for Her Quattro [™] Air / Quattro [™] Air for Her	Mirage Liberty [™]	Ultra Mirage [™]	Hospital Full Face
VPAP [™] Malibu	N/V	N/V	N/V	N/V	FULL FACE	FULL FACE	FULL FACE	FULL FACE	FULL FACE	FULL FACE
VPAP [™] S / VPAP [™] ST	N/V	N/V	N/V	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL	MIR FULL
VPAP [™] S (S9) / VPAP [™] ST (S9)	N/V	Full Face	N/V	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face
VPAP [™] ST-A (S9 iVAPS)	N/V	N/V	N/V	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face
VPAP [™] COPD (S9)	N/V	N/V	N/V	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face
S9 VPAP™ Tx	N/V	Full Face ²	N/V	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face
VPAP [™] Tx	N/V	N/V	N/V	N/V	FULL FACE	FULL FACE	FULL FACE	FULL FACE	FULL FACE	FULL FACE
Lumis [™] Tx	N/V	N/V	N/V	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face	Full Face
Stellar [™]	N/V	N/V	N/V	FULL FACE	FULL FACE	FULL FACE	FULL FACE	FULL FACE	FULL FACE	FULL FACE
VS Ultra [™] / VS Integra [™]	~	N/V	N/V	N/V	✓	~	N/V	~	~	~
VS III [™]	N/V	N/V	N/V	N/V	~	~	N/V	~	~	~
Astral [™]	FULL FACE	N/V	N/V	FULL FACE	N/V	FULL FACE	FULL FACE	FULL FACE	FULL FACE	FULL FACE

¹ SmartStop may not operate as described when using this mask with some CPAP or bi-level devices.

Nasal Masks 1

Device	AirFit [™] N30 for AirMini	AirFit [™] N30i ¹ / AirFit [™] N30	AirFit [™] N20 / AirFit [™] N20 for Her	AirTouch [™] N20 / AirTouch [™] N20 for Her	AirFit [™] N10 / AirFit [™] N10 for Her	Swift [™] FX Nano / Swift [™] FX Nano for Her	Mirage [™] FX / Mirage [™] FX for Her	Hospital Nasal
AirSense [™] 11	Χ	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal
AirMini™	~	Х	~	✓	Χ	Χ	Χ	Χ
AirStart [™] 10 / AirSense [™] 10	Χ	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal
AirCurve [™] 10 VAuto / AirCurve [™] 10 S / AirCurve [™] 10 ST	Х	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal
AirCurve [™] 10 ASV / AirCurve [™] 10 ASVAuto	Х	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal
AirCurve [™] 10 ST-A	Χ	N/V	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal
S9 AutoSet [™] / S9 AutoSet [™] for Her / S9 Elite [™]	Х	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal
S9 Escape [™] / S9 Escape [™] Auto	Χ	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal
C-Series Tango™	Χ	N/V	N/V	N/V	N/V	N/V	No setting required. Consult your clinician	before changing to a different mask type
S7 Elite [™]	Χ	N/V	N/V	N/V	N/V	N/V	STANDARD	STANDARD
S8 AutoSet Vantage [™]	Χ	N/V	SWIFT	SWIFT	SWIFT	SWIFT	STANDARD	STANDARD
S8 Compact [™]	Χ	N/V	N/V	N/V	No	setting required. Consult you	clinician before changing to a different	ent mask type.
S8 Elite [™] / S8 Elite [™] II	Χ	N/V	SWIFT	SWIFT	SWIFT	SWIFT	STANDARD	STANDARD
S8 Escape [™]	Χ	N/V	N/V	N/V	Require	es pressure calibration at 12 c	m H2O. Consult your clinician and che	eck the user guide.
S8 Escape [™] II	Χ	N/V		Requires	pressure calibration at 1	2 cm H2O. Consult your clinic	ian and check the user guide.	
S8 AutoSet [™] II	Χ	N/V	SWIFT	SWIFT	SWIFT	SWIFT	STANDARD	STANDARD
VPAP™ III	Χ	N/V	N/V	N/V	MIRAGE	MIRAGE	STANDARD	N/V
VPAP™ III ST-A	Х	N/V	N/V	N/V	MIRAGE	MIRAGE	STANDARD	N/V

² Excluding PAC and iVAPS modes. Only compatible for pressure range 4 - 25 cm H2O.

N/V = The mask has not been validated for this device.

 $[\]checkmark$ = Compatible X = Not Compatible

Device	AirFit [™] N30 for AirMini	AirFit [™] N30i ¹ / AirFit [™] N30	AirFit [™] N20 / AirFit [™] N20 for Her	AirTouch [™] N20 / AirTouch [™] N20 for Her	AirFit [™] N10 / AirFit [™] N10 for Her	Swift [™] FX Nano / Swift [™] FX Nano for Her	Mirage [™] FX / Mirage [™] FX for Her	Hospital Nasal
VPAP™ Adapt SV	Х	N/V	N/V	N/V	N/V	N/V	VISTA1	VISTA
VPAP™ Adapt (S9)	X	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal
VPAP™ Adapt (S9 ASV)	Х	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal
VPAP [™] Auto 25	Х	N/V	SWIFT	SWIFT	SWIFT	SWIFT	STANDARD	STANDARD
VPAP [™] Auto (S9)	Х	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal
VPAP [™] Malibu	Х	N/V	N/V	N/V	SWIFT	SWIFT	VISTA1	MERIDIAN
VPAP [™] S / VPAP [™] ST	Х	N/V	SWIFT	SWIFT	SWIFT	SWIFT	STANDARD	STANDARD
VPAP [™] S (S9) / VPAP [™] ST (S9)	Х	Pillows	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal
VPAP [™] ST-A (S9 iVAPS)	Х	N/V	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal
VPAP [™] COPD (S9)	Х	N/V	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal
S9 VPAP™ Tx	Х	Pillows2	Pillows	N/V	Pillows	Pillows	Nasal	Nasal
VPAP [™] Tx	Х	N/V	N/V	N/V	PILLOWS	PILLOWS	NASAL	NASAL
Lumis [™] Tx	Х	N/V	Pillows	Pillows	Pillows	Pillows	Nasal	Nasal
Stellar [™]	X	N/V	PILLOWS	PILLOWS	PILLOWS1	PILLOWS1	NASAL	NASAL
VS Ultra [™] / VS Integra [™]	X	N/V	N/V	N/V	N/V	N/V	N/V	N/V
VS III [™]	X	N/V	N/V	N/V	N/V	N/V	N/V	✓
Astral [™]	X	N/V	Pillows	Pillows	N/V	N/V	Nasal	Nasal

¹ SmartStop may not operate as described when using this mask with some CPAP or bi-level devices.

Nasal Masks 2

Device	Mirage Micro [™]	Ultra Mirage [™] II	Mirage Activa [™]	Mirage Activa [™] LT	Mirage [™] SoftGel	Mirage Vista [™]	Meridian	Mirage Kidsta [™]	Mirage Micro [™] for Kids	Pixi [™]
AirSense [™] 11	Х	X	X	Nasal	X	Х	Х	Χ	X	Х
AirMini™	Χ	Χ	Χ	Χ	Χ	Χ	Χ	Χ	Χ	Χ
AirStart [™] 10 / AirSense [™] 10	Nasal	Nasal	Nasal	Nasal	Nasal	Nasal	Nasal	Nasal	N/V	N/V
AirCurve [™] 10 VAuto / AirCurve [™] 10 S / AirCurve [™] 10 ST	Nasal	Nasal	Nasal	Nasal	Nasal	Nasal	Nasal	Nasal	N/V	N/V
AirCurve [™] 10 ASV / AirCurve [™] 10 ASVAuto	Nasal	Nasal	Nasal	Nasal	Nasal	Nasal	Nasal	Nasal	N/V	N/V
AirCurve [™] 10 ST-A	Nasal	Nasal	Nasal	Nasal	Nasal	Nasal	Nasal	Nasal	N/V	N/V
S9 AutoSet [™] / S9 AutoSet [™] for Her / S9 Elite [™]	Nasal	Nasal Ultra	Nasal	Nasal	Nasal	Nasal	Nasal	N/V	N/V	N/V
S9 Escape [™] / S9 Escape [™] Auto	Nasal	Nasal Ultra	Nasal	Nasal	Nasal	Nasal	Nasal	N/V	N/V	N/V
C-Series Tango [™]		No setting required.	Consult your clinician	before changing to a	different mask type.		N/V	N/V	N/V	N/V
S7 Elite [™]	STANDARD	ULTRA	STANDARD	STANDARD	STANDARD	STANDARD	N/V	N/V	N/V	N/V
S8 AutoSet Vantage [™]	STANDARD	ULTRA	ACTIVA	ACTIVA	ACTIVA	STANDARD	STANDARD	N/V	N/V	N/V
S8 Compact [™]		No settin	g required. Consult yo	ur clinician before ch	anging to a different ma	ask type.		N/V	N/V	N/V
S8 Elite [™] / S8 Elite [™] II	STANDARD	ULTRA	ACTIVA	ACTIVA	ACTIVA	STANDARD	STANDARD	N/V	N/V	N/V
S8 Escape [™]		Requires pres	ssure calibration at 12	cm H2O. Consult you	r clinician and check th	e user guide.		N/V	N/V	N/V

² Excluding PAC and iVAPS modes. Only compatible for pressure range 4 - 25 cm H2O.

N/V = The mask has not been validated for this device.

 $[\]checkmark$ = Compatible X = Not Compatible

Device	Mirage Micro [™]	Ultra Mirage [™] II	Mirage Activa [™]	Mirage Activa [™] LT	Mirage [™] SoftGel	Mirage Vista [™]	Meridian	Mirage Kidsta [™]	Mirage Micro [™] for Kids	Pixi™
S8 Escape [™] II		Requires pres	ssure calibration at 12	cm H2O. Consult you	r clinician and check th	e user guide.		N/V	N/V	N/V
S8 AutoSet [™] II	STANDARD	ULTRA	ACTIVA	ACTIVA	ACTIVA	STANDARD	STANDARD	N/V	N/V	N/V
VPAP™ III	STANDARD	ULTRA	STANDARD	STANDARD	STANDARD	STANDARD	N/V	N/V	N/V	N/V
VPAP™ III ST-A	STANDARD	ULTRA	STANDARD	STANDARD	STANDARD	STANDARD	N/V	STANDARD	STANDARD	N/V
VPAP™ Adapt SV	VISTA1	ULTRA	ACTIVA	ACTIVA	ACTIVA	VISTA	VISTA	N/V	N/V	N/V
VPAP [™] Adapt (S9)	Nasal	Nasal Ultra	Nasal	Nasal	Nasal	Nasal	Nasal	N/V	N/V	N/V
VPAP™ Adapt (S9 ASV)	Nasal	Nasal Ultra	Nasal	Nasal	Nasal	Nasal	Nasal	N/V	N/V	N/V
VPAP [™] Auto 25	STANDARD	ULTRA	ACTIVA	ACTIVA	ACTIVA	STANDARD	STANDARD	N/V	N/V	N/V
VPAP [™] Auto (S9)	Nasal	Nasal Ultra	Nasal	Nasal	Nasal	Nasal	Nasal	N/V	N/V	N/V
VPAP [™] Malibu	MERIDIAN	ULTRA	N/V	N/V	N/V	VISTA	MERIDIAN	N/V	N/V	N/V
VPAP™S / VPAP™ST	STANDARD	ULTRA	ACTIVA	ACTIVA	ACTIVA	STANDARD	STANDARD	N/V	N/V	N/V
VPAP [™] S (S9) / VPAP [™] ST (S9)	Nasal	Nasal Ultra	Nasal	Nasal	Nasal	Nasal	Nasal	N/V	N/V	N/V
VPAP™ ST-A (S9 iVAPS)	Nasal	Nasal Ultra	Nasal	Nasal	Nasal	Nasal	Nasal	N/V	N/V	Pediatric
VPAP [™] COPD (S9)	Nasal	Nasal Ultra	Nasal	Nasal	Nasal	Nasal	Nasal	N/V	N/V	Pediatric
S9 VPAP™ Tx	Nasal	Nasal Ultra	Nasal	Nasal	Nasal	Nasal	Nasal	N/V	N/V	Pediatric
VPAP [™] Tx	NASAL	NASAL-ULTRA	NASAL	NASAL	NASAL	NASAL	N/V	NASAL	N/V	N/V
Lumis [™] Tx	Nasal	Nasal Ultra	Nasal	Nasal	Nasal	Nasal	Nasal	N/V	N/V	Pediatric
Stellar [™]	NASAL	ULTRA	NASAL	NASAL	NASAL	NASAL	N/V	NASAL	N/V	Pediatric
VS Ultra [™] / VS Integra [™]	N/V	N/V	N/V	N/V	N/V	N/V	N/V	N/V	N/V	N/V
VS III [™]	✓	~	~	~	N/V	~	N/V	~	N/V	N/V
Astral [™]	Nasal	Nasal Ultra	N/V	Nasal	Nasal	Nasal	N/V	Nasal	Nasal	Pediatric

N/V = The mask has not been validated for this device.

 \checkmark = Compatible X = Not Compatible

Nasal Pillows

Device	AirFit [™] P30i ¹	AirFit [™] P10 For AirMini [™]	AirFit [™] P10 / AirFit [™] P10 for Her	Swift [™] FX / Swift [™] FX for Her	Swift [™] FX Bella / Swift [™] FX Bella Gray	Swift [™] LT / Swift [™] LT for Her	Mirage Swift [™] II
AirSense [™] 11	Pillows	Χ	Pillows2	Pillows	Pillows	Pillows	Pillows
AirMini™	Χ		Х	X	Χ	Х	Х
AirStart [™] 10 / AirSense [™] 10	Pillows	Χ	Pillows2	Pillows	Pillows	Pillows	Pillows
AirCurve [™] 10 VAuto / AirCurve [™] 10 S / AirCurve [™] 10 ST	Pillows	X	Pillows2	Pillows	Pillows	Pillows	Pillows
AirCurve [™] 10 ASV / AirCurve [™] 10 ASVAuto	Pillows	X	Pillows2	Pillows	Pillows	Pillows	Pillows
AirCurve [™] 10 ST-A	N/V	Χ	Pillows2	Pillows	Pillows	Pillows	Pillows
S9 AutoSet [™] / S9 AutoSet [™] for Her / S9 Elite [™]	Pillows	X	Pillows	Pillows	Pillows	Pillows	Pillows
S9 Escape [™] / S9 Escape [™] Auto	Pillows	Χ	Pillows	Pillows	Pillows	Pillows	Pillows
C-Series Tango [™]	N/V	Χ		No setting required. Cons	sult your clinician before changing	to a different mask type.	
S7 Elite [™]	N/V	Χ	MIRAGE	MIRAGE	MIRAGE	MIRAGE	MIRAGE
S8 AutoSet Vantage [™]	N/V	Х	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT

Device	AirFit [™] P30i ¹	AirFit [™] P10 For AirMini [™]	AirFit [™] P10 / AirFit [™] P10 for Her	Swift [™] FX / Swift [™] FX for Her	Swift [™] FX Bella / Swift [™] FX Bella Gray	Swift [™] LT / Swift [™] LT for Her	Mirage Swift [™] II			
S8 Compact [™]	N/V	Х		No setting required. Consult your clinician before changing to a different mask type.						
S8 Elite [™] / S8 Elite [™] II	N/V	Х	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT			
S8 Escape [™]	N/V	Χ		Requires pressure calibration	n at 12 cm H2O. Consult your clinic	cian and check the user guide.				
S8 Escape [™] II	N/V	Х		Requires pressure calibration	n at 12 cm H2O. Consult your clinic	cian and check the user guide.				
S8 AutoSet [™] II	N/V	Х	SWIFT	SWIFT	SWIFT	SWIFT	SWIFT			
VPAP™ III	N/V	Х	MIRAGE ²	MIRAGE	MIRAGE	MIRAGE	MIRAGE			
VPAP™ III ST-A	N/V	Х	MIRAGE ²	MIRAGE	MIRAGE	MIRAGE	MIRAGE			
VPAP [™] Adapt SV	N/V	Х	N/V	N/V	N/V	N/V	N/V			
VPAP [™] Adapt (S9)	Pillows	Χ	Pillows ²	Pillows	Pillows	Pillows	Pillows			
VPAP [™] Adapt (S9 ASV)	Pillows	Χ	Pillows ²	Pillows	Pillows	Pillows	Pillows			
VPAP [™] Auto 25	N/V	Х	SWIFT ²	SWIFT	SWIFT	SWIFT	SWIFT			
VPAP [™] Auto (S9)	Pillows	Χ	Pillows ²	Pillows	Pillows	Pillows	Pillows			
VPAP™ Malibu	N/V	Χ	SWIFT ²	SWIFT	SWIFT	SWIFT	SWIFT			
$VPAP^{M}S/VPAP^{M}ST$	N/V	Х	SWIFT ²	SWIFT	SWIFT	SWIFT	SWIFT			
$VPAP^{M} S (S9) / VPAP^{M} ST (S9)$	Pillows	Х	Pillows ²	Pillows	Pillows	Pillows	Pillows			
VPAP™ ST-A (S9 iVAPS)	N/V	Х	Pillows ²	Pillows	Pillows	Pillows	Pillows			
VPAP™ COPD (S9)	N/V	Χ	Pillows ²	Pillows	Pillows	Pillows	Pillows			
S9 VPAP™ Tx	Pillows3	Х	Pillows ²	Pillows	Pillows	Pillows	Pillows			
VPAP™ Tx	N/V	Χ	PILLOWS ²	PILLOWS	PILLOWS	PILLOWS	PILLOWS			
Lumis [™] Tx	N/V	Х	Pillows ²	Pillows	Pillows	Pillows	Pillows			
Stellar [™]	N/V	Χ	PILLOWS ^{1,2}	PILLOWS ¹	PILLOWS ¹	PILLOWS ¹	PILLOWS ¹			
VS Ultra [™] / VS Integra [™]	N/V	Χ	N/V	N/V	N/V	N/V	N/V			
VS III [™]	N/V	Χ	N/V	N/V	N/V	N/V	N/V			
Astral [™]	N/V	Х	N/V	Pillows	Pillows	Pillows	Pillows			

¹ SmartStop may not operate as described when using this mask with some CPAP or bi-level devices.

N/V = The mask has not been validated for this device.

 \checkmark = Compatible X = Not Compatible



² Accuracy of tidal volume/minute ventilation may reduce using the AirFit P10 mask in bilevel modes (S, T, ST), particularly at higher pressure supports.

³ Excluding PAC and iVAPS modes. Only compatible for pressure range 4 - 25 cm H2O.